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Dmuchana miska do mycia, ADVANCE

Model: PM200AM

Wstęp

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Model PM200AM

Introduction

Thank you for choosing our product. We are confident that it will meet your expectations. We encourage you to read this document before using the product for the first time. These instructions are intended to provide the necessary information on the safe use of the medical device. The manual provides detailed guidance on the device's functions and capabilities. Retaining the document allows this information to be easily found in the future, facilitating correct use.

Description of the medical device

A medical device designed for washing the head of people lying down with a gravity shower, who cannot use a traditional bathroom due to a medical condition. The bowl is made of a durable, leak-proof material (usually PVC or vinyl), has two air chambers, a special neck cut-out for comfort and protection against splashing, and a drain for easy removal of water after washing. It is designed for the daily care of the elderly, disabled, with limited mobility or chronically lying down, making it significantly easier for carers to maintain the patient's hygiene without having to move them from bed. Thanks to its pneumatic design, the bowl is lightweight, easy to store and transport, and its use minimises the risk of bedding getting wet and increases the comfort of both patient and carer's comfort.

Intended use

The product is designed for the head and hair hygiene of patients immobilised in bed, including the elderly, the disabled, those undergoing surgery or those suffering from chronic illnesses that prevent independent movement to the bathroom. The product enables hair to be washed safely and comfortably without having to change the patient's position, which is particularly important in long-term care, hospice care and geriatric hospitals and rehabilitation wards.

Contraindications

Allergy to device materials, especially PVC or vinyl, which can cause contact allergic reactions

in sensitive patients. In addition, the device should not be used in patients with open wounds or skin lesions in the neck and head area, where there is a risk of infection or wounds being aggravated by contact with water and bowl material.

Description of product construction

The inflatable wash bowl is an ergonomically designed medical device in white, made of soft, inflatable PVC material in a distinctive horseshoe shape with a deep neckline. The design consists of two main air chambers which, when inflated, form a stable and comfortable base, equipped with a special drain hole located at the lowest point of the bowl for controlled drainage of used water. The surface of the bowl is smooth, easy to clean and resistant to disinfectants, and thanks to its pneumatic design, it takes up little space when deflated, making it easy to store and transport. The gravity shower version is additionally equipped with a flexible water reservoir with a hanging bracket, a conduit and a mini shower head, which makes it possible to conveniently rinse the lying patient's hair with free-flowing water.

Installation / first use

Preparing the bowl

- Inflate the bowl using your mouth or the supplied air pump.
- Make sure the bowl is properly inflated and stable.

Preparation of the drainage system

- Close the drain valve in the bowl.
- Attach the drain hose to the drain valve.
- Place the other end of the drain hose in the waste water container, which must be below the level of the bowl.

Preparing the gravity shower

- Hang the clean water tank at a suitable height, such as on a stand or hook above the level of the bowl.
- Connect the conduit with the mini shower head to the water tank.
- Check the operation of the valve on the shower head - water should flow freely under the force of gravity when it is opened.

Bowl positioning

- Place the inflatable bowl under the patient's head using the special neck cut-out.
- Make sure the bowl is firmly placed on a hard surface.
- Check that the drain hose is not pinched or kinked.

The user or their carer should consult a healthcare professional if they are unsure whether the medical device is appropriate for the patient in question.

Any use not in accordance with the intended use is considered inappropriate and DANGEROUS. The manufacturer and/or distributor shall not be liable for damage caused by improper, incorrect and/or unreasonable use.

Safety and usage rules

- Before filling with water, first pour in cold, then warm water; the water temperature should not exceed 40°C.
- Check the bowl for damage, cracks or leaks before each use.
- Ensure that all valves and connections are tight and stable.
- Do not use the bowl near sharp objects that may damage it.
- Place the bowl under the patient's head so that the neckline provides comfort and prevents water leakage.
- Check that the drain hose is properly secured and directed to the waste water container below the level of the bowl.
- Only use the bowl on a flat, stable surface.
- Always monitor the temperature of the water - it should be warm but not hot.
- Do not leave the patient unattended during washing.
- Ensure that the drain hose is not kinked or pinched to allow water to drain freely.
- When you have finished washing, empty the bowl of water, clean and disinfect it with suitable products, and allow it to dry before folding and storing.

Maintenance

1. Regular cleaning of the product should be carried out with water, a sponge and mild detergents intended for plastics, such as soap or a mild cleaning agent approved by the National Institute of Hygiene. It is important to dry the product thoroughly after each use.
2. Avoid the use of: strong detergents, abrasives, bleach, petrol and all solvents, as these may cause damage to the product, visible defects and lead to a loss of warranty.
3. After each use, the device should be disinfected using disinfectants designed for disinfection of surfaces in contact with the skin, according to the manufacturer's recommendations, and then thoroughly rinsed and dried to prevent the growth of bacteria or fungi.
4. Store the product in a dry, ventilated place, away from heat sources and corrosive substances; do not fold or store the product in a damp state to avoid mould formation and material damage.
5. Regularly check the condition of the product, paying attention to the tightness of the air chambers and the absence of mechanical damage; if minor leaks are detected, the repair kit supplied with the product can be used.

Storage and transport conditions

The product should be stored in a dry and well-ventilated place, protected from moisture and extreme temperatures to prevent damage to the materials. Contact with chemicals that may affect its durability and performance should be avoided. During transport, the product should be secured in such a way as to prevent displacement and damage. Use appropriate protective packaging.

Disposal

Disposal of the product must be carried out in accordance with the applicable local environmental regulations.

Product life

The expected product life is up to 5 years if used in accordance with the information in the instructions for use. The product can be used after this period if it is undamaged and meets the specified requirements for safe use.

Warranty

- The product is covered by a 24-month warranty from the date of purchase, against any defects in material or workmanship.
- The warranty covers the free replacement and/or repair of parts found to be manufacturing defects.
- The warranty does not cover parts subject to normal wear and tear.
- The warranty is no longer valid if the product has been tampered with/modified, if the defect is the result of misuse or if the damage cannot be attributed to the manufacturer.

NOTE

Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is residing.

Product content

The product is made of high-quality, durable PVC (polyvinyl chloride), with a thickness typically in the range 0.3-0.5 mm, which ensures flexibility, resistance to mechanical damage and ease of hygiene.

Technical data

- External dimensions: approx. 43 x 36 x 39 cm
- Gravity shower set
- Weight: approx. 0.5 kg

Explanation of symbols

	Conformity of the device with EU Regulation 2017/745 on medical devices		Medical device
	Take a look at the instructions for use		Batch/series number
	Model		Warning
	Date of production		Medical device manufacturer

	Unique identification code of the device		Keep out of sunlight
	Protect from moisture		

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ACTION S.A.